

# MEDCASTER PT/PS/PC





# MEDCASTER PRIMACY PT/PS/PC SERIES

Prime solutions for your hospital beds, stretchers and diagnostic equipment










Ask about our free hospital guide and cross reference

# MEDCASTER PRIMACY PT/PS/PC SERIES

Prime solutions for your hospital beds, stretchers and diagnostic equipment

- QUIET
- CLEAN
- Ѕмоотн
- MAINTENANCE FREE



### Models

#### **PT Twin Wheel Models**

- Wheel diameter models include:
- 100mm (125kg/275lbs capacity each)
- 125mm (130kg/286lbs capacity each)
- 150mm (150kg/330lbs capacity each)

#### **PS Single Wheel Models**

• 125mm, 150mm and 200mm wheel diameter models


#### **PC Single Wheel Models**

- 100mm, 125mm, 150mm and 200mm wheel diameter models
- Optional conductive models
- Optional colors with quantity orders

## **Applications**

- Care beds
- Hospital beds
- Intensive care beds
- Medical appliances
- OR beds
- OR tables
- Stretchers
- Gurneys
- Medical recliners
- Mobile ultrasound equipment

# **Mounting Options**



- Stem diameter 27.8mmTwo locking hole options
- availableLocking hole in line or across direction of travel
- Four fixing holes
- Snap-on adapter sleeves available for different tube dimensions
- Standard cam (41°/52°)

#### CL02

- Stem diameter 32mm
- Two locking hole options available
- Locking hole in line or across direction of travel
- Four fixing holes
- Standard cam 30°

# Pedal/Lever Options (7/16' hex)



Transmission Levers



866-462-9700

WWW.MEDCASTER.COM
New Size! 399 Himbus

Available in 150mm – New 125mm wheel diameter sizes coming December 2011.



Availa New comin

